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B'MAGIE™ intelligent circulation booster

Quick Guide

Please read the User Manual for details information & follow your doctor's instruction.

1

Put the B'MAGIE™ in front of your seat, and connect to the socket nearby.

2

Seat on a chair, with the feet falling down, then step on the pad of B'MAGIE™ gently.

3

Press the “” button on panel or remote control, to start the session.

Precautions and Contra-Indication:

- Don't stand on the device, enjoy your session on chair
- Don't splash any liquid on the device, or use it in moisture environment
- Pregnancy woman is forbidden
- Don't use it on the foot with infection or wound
- Adult only

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After service: Shaanxi Bachmen Health Product Ltd.

Address: B-1 Room, 35 Tangya Road, Xiaan Hi-Tech Zone, China

Tel: + 86 400 600 9971

4

The default time “30”minutes shows on the panel, adjust the timer with “+” or “-” (step by 1 minuter).

5

Same time, the intension shows “0” beside the sign of , you can adjust it by “+” or “-” between 0~99. (Longer press the button could adjust it by step of 10, short press step of 1) Modify it until you could tolerance the intension.

6

When the countdown timer go to “0”, B'MAGIE™ stopped automatically. You can interrupt the session whenever you press the “” button on the panel or remote control.

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Authorization: Bachmen Gesundheit Produkt GmbH